

KiRA

Kinematic Rapid Assessment



THE ULTIMATE
ACL SUPERVISOR

ORTHOKEY introduces a new generation of smart instruments for orthopaedics.

A RELIABLE COMPANION


KiRA quantifies the dynamic laxity of the knee during Pivot shift test and Lachman test. Easy to use, non-invasive, small sized and portable, KiRA supports an accurate diagnosis in every context.

KiRA is the first and exclusive digital medical device CE certified to measure knee laxity.


KiRA is composed by a free App available on Apple Store and by a sensing device that communicates via Bluetooth® technology with iPad 3 or superior, iPad mini and iPhone 5 or superior.

WITH KiRA YOU CAN


- Detect and assess the grade of suspected ACL lesion
- Analyze the functional recovery during rehabilitation
- Patient's history analysis for an effective follow up evaluation



1



2



3

Advantages for the professional



Fast, Effective, Immediate

- Provides an exact assessment and immediate visualization of the knee stability grade level
- Crucial to define the best clinical approach in case of partial lesion

Advantages for the patient



Capture your knee health condition

- Helps you to comprehend the diagnosis
- Tracks your improvements during rehabilitation

Advantages for both



Greater security in diagnosis!

- Reduces false positive risks through standardized diagnostic process
- Improves the level of confidence of doctor and patient



This document is the confidential property of ORTHOKEY. No authorised reproduction and/or use of any portion is allowed without prior written approval. Pictures are for demonstration purposes only. Orthokey reserves the right to make changes at any time to the interfaces without prior notice.

www.smartkira.com

KiRA.BREN.1509